

# IPPAI

version 0.1.0b

**2020, CAMI (Computer Assisted Medical Interventions), DKFZ, Heidelberg**

April 01, 2020



# Contents

<b>Welcome to IPPAI's documentation!</b>	<b>1</b>
<b>README</b>	<b>1</b>
Internal Install Instructions	1
Building the documentation	1
Overview	1
Simulating photoacoustic images	1
<b>Examples</b>	<b>1</b>
Performing a simple optical forward simulation	2
Reading the HDF5 simulation output	3
Defining custom tissue structures and properties	4
<b>Class references</b>	<b>5</b>
<b>Index</b>	<b>7</b>
<b>Python Module Index</b>	<b>9</b>



# Welcome to IPPAI's documentation!



## README

The Image Processing for Photoacoustic Imaging (IPPAI) toolkit.

## Internal Install Instructions

These install instructions are made under the assumption that you have access to the phabricator ippai project. When you are reading these instructions there is a 99% chance that is the case (or someone send these instructions to you).

So, for the 1% of you: Please also follow steps 1 - 3:

1. `git clone https://phabricator.mtk.org/source/ippai.git`
2. `git checkout master`
3. `git pull`

Now open a python instance in the 'ippai' folder that you have just downloaded. Make sure that you have your preferred virtual environment activated

1. `cd ippai`
2. `python -m setup.py build install`
3. Test if the installation worked by using `python` followed by `import ippai` then `exit()`

If no error messages arise, you are now setup to use ippai in your project.

## Building the documentation

When the installation went fine and you want to make sure that you have the latest documentation you should do the following steps in a command line:

1. Navigate to the ippai source directory (same level where the setup.py is in)
2. Execute the command `sphinx-build -b pdf -a documentation/src documentation`
3. Find the PDF file in `documentation/ippai_documentation.pdf`

## Overview

The main use case for the ippai framework is the simulation of photoacoustic images. However, it can also be used for image processing.

## Simulating photoacoustic images

A basic example on how to use ippai in you project to run an optical forward simulation is given in the `samples/minimal_optical_simulation.py` file.

## Examples

## Performing a simple optical forward simulation

The file can be found in `samples/minimal_optical_simulation.py`:

```
from ippai.utils import Tags

from ippai.simulate.simulation import simulate
from ippai.simulate.structures import create_epidermis_layer
from ippai.simulate.structures import create_muscle_background
from ippai.simulate.structures import create_vessel_tube

import numpy as np

# TODO change these paths to the desired executable and save folder
SAVE_PATH = "path/to/save/file"
MCX_BINARY_PATH = "path/to/mcx/binary"

VOLUME_WIDTH_IN_MM = 10
VOLUME_HEIGHT_IN_MM = 10
SPACING = 0.25
RANDOM_SEED = 4711

def create_example_tissue():
    """
    This is a very simple example script of how to create a tissue definition.
    It contains a muscular background, an epidermis layer on top of the muscles
    and a blood vessel.
    """
    tissue_dict = dict()
    tissue_dict["background"] = create_muscle_background()
    tissue_dict["epidermis"] = create_epidermis_layer()
    tissue_dict["vessel"] = create_vessel_tube(x_min=0, x_max=VOLUME_WIDTH_IN_MM,
                                              z_min=0, z_max=VOLUME_HEIGHT_IN_MM,
                                              r_min=1, r_max=3)

    return tissue_dict

# Seed the numpy random configuration prior to creating the settings file in
# order to ensure that the same volume
# is generated with the same random seed every time.

np.random.seed(RANDOM_SEED)

settings = {
    # These parameters set the general properties of the simulated volume
    Tags.RANDOM_SEED: RANDOM_SEED,
    Tags.VOLUME_NAME: "MyVolumeName_" + str(RANDOM_SEED),
    Tags.SIMULATION_PATH: SAVE_PATH,
    Tags.SPACING_MM: SPACING,
    Tags.DIM_VOLUME_Z_MM: VOLUME_HEIGHT_IN_MM,
    Tags.DIM_VOLUME_X_MM: VOLUME_WIDTH_IN_MM,
    Tags.DIM_VOLUME_Y_MM: VOLUME_WIDTH_IN_MM,
    Tags.AIR_LAYER_HEIGHT_MM: 0,
    Tags.GELPAD_LAYER_HEIGHT_MM: 0,

    # The following parameters set the optical forward model
    Tags.RUN_OPTICAL_MODEL: True,
    Tags.WAVELENGTHS: np.arange(700, 951, 10),
    Tags.OPTICAL_MODEL_NUMBER_PHOTONS: 1e7,
    Tags.OPTICAL_MODEL_BINARY_PATH: MCX_BINARY_PATH,
    Tags.OPTICAL_MODEL: Tags.MODEL_MCX,
```

```

Tags.ILLUMINATION_TYPE: Tags.ILLUMINATION_TYPE_PENCIL,
Tags.LASER_PULSE_ENERGY_IN_MILLIJOULE: 50,

# The following parameters tell the script that we do not want any extra
# modelling steps
Tags.RUN_ACOUSTIC_MODEL: False,
Tags.APPLY_NOISE_MODEL: False,
Tags.PERFORM_IMAGE_RECONSTRUCTION: False,
Tags.SIMULATION_EXTRACT_FIELD_OF_VIEW: False,

# Add the structures to be simulated to the tissue
Tags.STRUCTURES: create_example_tissue()
}
print("Simulating ", RANDOM_SEED)
simulate(settings)
# TODO settings[Tags.IPPAI_OUTPUT_PATH]
print("Simulating ", RANDOM_SEED, "[Done]")

```

## Reading the HDF5 simulation output

The file can be found in `samples/access_saved_PA1_data.py`:

```

from ippai.io_handling import load_hdf5, save_hdf5
import matplotlib.pyplot as plt
import numpy as np
from ippai.simulate import SaveFilePaths

PATH = "/home/janek/test/Vessels_10005/ippai_output.hdf5"
WAVELENGTH = 800 # currently only 800 and 900 are simulated as well

file = load_hdf5(PATH)

print(file['simulations'].keys())

fluence = (file['simulations']['original_data']['optical_forward_model_output']
           [str(WAVELENGTH)]['fluence'])
initial_pressure = (file['simulations']['original_data']
                    ['optical_forward_model_output']
                    [str(WAVELENGTH)]['initial_pressure'])
absorption = (file['simulations']['original_data']['simulation_properties']
              [str(WAVELENGTH)]['mua'])

shape = np.shape(fluence)

if len(shape) > 2:
    plt.figure()
    plt.subplot(231)
    plt.imshow(np.log10(fluence[int(shape[0]/2), :, :]))
    plt.subplot(232)
    plt.imshow(np.log10(absorption[int(shape[0]/2), :, :]))
    plt.subplot(233)
    plt.imshow(np.log10(initial_pressure))
    plt.subplot(234)
    plt.imshow(np.log10(fluence[:, int(shape[1]/2), :]))
    plt.subplot(235)
    plt.imshow(np.log10(absorption[:, int(shape[1]/2), :]))
    plt.subplot(236)
    plt.imshow(np.log10(initial_pressure))
    plt.show()
else:

```

```

plt.figure()
plt.subplot(131)
plt.imshow(np.log10(fluence[1:129, -65:-1]))
plt.subplot(132)
plt.imshow(np.log10(absorption[1:129, -65:-1]))
plt.subplot(133)
plt.imshow(np.log10(initial_pressure[1:129, -65:-1]))
plt.show()

save_hdf5(file, _PATH)

```

## Defining custom tissue structures and properties

The file can be found in `samples/create_custom_tissues.py`:

```

from ippai.utils import TissueSettingsGenerator
from ippai.utils import CHROMOPHORE_LIBRARY
from ippai.utils import Chromophore
from ippai.utils import AbsorptionSpectrum
import numpy as np

def create_custom_absorber():
    wavelengths = np.linspace(200, 1500, 100)
    absorber = AbsorptionSpectrum(spectrum_name="random absorber",
                                wavelengths=wavelengths,
                                absorption_per_centimeter=np.random.random(
                                    np.shape(wavelengths)))

    return absorber

def create_custom_chromophore(volume_fraction: float = 1.0):
    chromophore = Chromophore(
        spectrum=create_custom_absorber(),
        volume_fraction=volume_fraction,
        musp500=40.0,
        b_mie=1.1,
        f_ray=0.9,
        anisotropy=0.9
    )
    return chromophore

def create_custom_tissue_type():
    # First create an instance of a TissueSettingsGenerator
    tissue_settings_generator = TissueSettingsGenerator()

    water_volume_fraction = 0.4
    bvf = 0.5
    oxy = 0.4

    # Then append chromophores that you want
    tissue_settings_generator.append(key="oxyhemoglobin", value=
        CHROMOPHORE_LIBRARY.oxyhemoglobin(oxy*bvf))
    tissue_settings_generator.append(key="deoxyhemoglobin", value=
        CHROMOPHORE_LIBRARY.deoxyhemoglobin(oxy * bvf))
    tissue_settings_generator.append(key="water", value=
        CHROMOPHORE_LIBRARY.water(water_volume_fraction))
    tissue_settings_generator.append(key="custom", value=

```



```

        create_custom_chromophore(0.1))

return tissue_settings_generator.get_settings()

```

## Class references

### `class ippai.utils.MorphologicalTissueProperties`

This class contains a listing of morphological tissue parameters as reported in literature. The listing is not the result of a meta analysis, but rather uses the best fitting paper at the time of implementation. Each of the fields is annotated with a literature reference or a descriptions of how the particular values were derived for tissue modelling.

### `class ippai.utils.OpticalTissueProperties`

This class contains a listing of optical tissue parameters as reported in literature. The listing is not the result of a meta analysis, but rather uses the best fitting paper at the time of implementation. Each of the fields is annotated with a literature reference or a descriptions of how the particular values were derived for tissue modelling.

### `class ippai.utils.StandardProperties`

This class contains a listing of default parameters that can be used. These values are sensible default values but are generally not backed up by proper scientific references, or are rather specific for internal use cases.

`class ippai.utils.AbsorptionSpectrum (spectrum_name: str, wavelengths: numpy.ndarray, absorption_per_centimeter: numpy.ndarray)`

An instance of this class represents the absorption spectrum over wavelength for a particular

`get_absorption_for_wavelength (wavelength: float) → float`

**Parameters:** **wavelength** – the wavelength to retrieve a optical absorption value for [cm<sup>-1</sup>].

**Returns:** the best matching linearly interpolated absorption value for the given wavelength.

`get_absorption_over_wavelength ()`

**Returns:** numpy array with the available wavelengths and the corresponding absorption properties

### `class ippai.utils.AbsorptionSpectrumLibrary`

`ippai.utils.view_absorption_spectra (save_path=None)`

Opens a matplotlib plot and visualizes the available absorption spectra.

**Parameters:** **save\_path** – If not None, then the figure will be saved as a png file to the destination.

`class ippai.utils.Chromophore (spectrum: ippai.utils.libraries.spectra_library.AbsorptionSpectrum, volume_fraction: float, musp500: float, f_ray: float, b_mie: float, anisotropy: float)`

### `class ippai.utils.ChromophoreLibrary`

### `class ippai.utils.TissueLibrary`

TODO

`blood_arterial ()`

**Returns:** a settings dictionary containing all min and max parameters fitting for full blood.

`blood_generic ()`

**Returns:** a settings dictionary containing all min and max parameters fitting for full blood.

`blood_venous ()`

**Returns:** a settings dictionary containing all min and max parameters fitting for full blood.

**bone ()**

**Returns:** a settings dictionary containing all min and max parameters fitting for full blood.

**constant (mua, mus, g)**  
TODO

**dermis (background\_oxy=0.5)**

**Returns:** a settings dictionary containing all min and max parameters fitting for dermis tissue.

**epidermis ()**

**Returns:** a settings dictionary containing all min and max parameters fitting for epidermis tissue.

**get\_blood\_volume\_fractions (total\_blood\_volume\_fraction, oxygenation)**  
TODO

**muscle (background\_oxy=0.5)**

**Returns:** a settings dictionary containing all min and max parameters fitting for generic background tissue.

**subcutaneous\_fat (background\_oxy=0.5)**

**Returns:** a settings dictionary containing all min and max parameters fitting for subcutaneous fat tissue.

**class ippai.utils.TissueSettingsGenerator**  
TODO

# Index

## A

AbsorptionSpectrum (class in ippai.utils)  
AbsorptionSpectrumLibrary (class in ippai.utils)

## B

blood\_arterial() (ippai.utils.TissueLibrary method)  
blood\_generic() (ippai.utils.TissueLibrary method)  
blood\_venous() (ippai.utils.TissueLibrary method)  
bone() (ippai.utils.TissueLibrary method)

## C

Chromophore (class in ippai.utils)  
ChromophoreLibrary (class in ippai.utils)  
constant() (ippai.utils.TissueLibrary method)

## D

dermis() (ippai.utils.TissueLibrary method)

## E

epidermis() (ippai.utils.TissueLibrary method)

## G

get\_absorption\_for\_wavelength()  
(ippai.utils.AbsorptionSpectrum method)  
get\_absorption\_over\_wavelength()  
(ippai.utils.AbsorptionSpectrum method)  
get\_blood\_volume\_fractions() (ippai.utils.TissueLibrary method)

## I

ippai.utils (module)

## M

MorphologicalTissueProperties (class in ippai.utils)  
muscle() (ippai.utils.TissueLibrary method)

## O

OpticalTissueProperties (class in ippai.utils)

## S

StandardProperties (class in ippai.utils)  
subcutaneous\_fat() (ippai.utils.TissueLibrary method)

## T

TissueLibrary (class in ippai.utils)  
TissueSettingsGenerator (class in ippai.utils)

## V

view\_absorption\_spectra() (in module ippai.utils)



# Python Module Index

## *i*

ippai

[ippai.utils](#)