

IPPAI

version 0.1.0

CAMI (Computer Assisted Medical Interventions), DKFZ, Heidelberg

April 01, 2020

Contents

Welcome to IPPAI's documentation!	1
README	1
Internal Install Instructions	1
Building the documentation	1
Overview	1
Simulating photoacoustic images	1
Examples	1
Performing a simple optical forward simulation	2
Reading the HDF5 simulation output	3
Defining custom tissue structures and properties	4
Class references	5
Module: utils	5
Module: io_handling	6
Module: simulate	6
Conventions	8
Structures	8
Forward models	8
Index	11
Python Module Index	13

Welcome to IPPAI's documentation!



README

The Image Processing for Photoacoustic Imaging (IPPAI) toolkit.

Internal Install Instructions

These install instructions are made under the assumption that you have access to the phabricator ippai project. When you are reading these instructions there is a 99% chance that is the case (or someone send these instructions to you).

So, for the 1% of you: Please also follow steps 1 - 3:

1. `git clone https://phabricator.mtk.org/source/ippai.git`
2. `git checkout master`
3. `git pull`

Now open a python instance in the 'ippai' folder that you have just downloaded. Make sure that you have your preferred virtual environment activated

1. `cd ippai`
2. `python -m setup.py build install`
3. Test if the installation worked by using `python` followed by `import ippai` then `exit()`

If no error messages arise, you are now setup to use ippai in your project.

Building the documentation

When the installation went fine and you want to make sure that you have the latest documentation you should do the following steps in a command line:

1. Navigate to the ippai source directory (same level where the setup.py is in)
2. Execute the command `sphinx-build -b pdf -a documentation/src documentation`
3. Find the PDF file in `documentation/ippai_documentation.pdf`

Overview

The main use case for the ippai framework is the simulation of photoacoustic images. However, it can also be used for image processing.

Simulating photoacoustic images

A basic example on how to use ippai in you project to run an optical forward simulation is given in the `samples/minimal_optical_simulation.py` file.

Examples

Performing a simple optical forward simulation

The file can be found in `samples/minimal_optical_simulation.py`:

```
from ippai.utils import Tags

from ippai.simulate.simulation import simulate
from ippai.simulate.structures import create_epidermis_layer
from ippai.simulate.structures import create_muscle_background
from ippai.simulate.structures import create_vessel_tube

import numpy as np

# TODO change these paths to the desired executable and save folder
SAVE_PATH = "path/to/save/file"
MCX_BINARY_PATH = "path/to/mcx/binary"

VOLUME_WIDTH_IN_MM = 10
VOLUME_HEIGHT_IN_MM = 10
SPACING = 0.25
RANDOM_SEED = 4711

def create_example_tissue():
    """
    This is a very simple example script of how to create a tissue definition.
    It contains a muscular background, an epidermis layer on top of the muscles
    and a blood vessel.
    """
    tissue_dict = dict()
    tissue_dict["background"] = create_muscle_background()
    tissue_dict["epidermis"] = create_epidermis_layer()
    tissue_dict["vessel"] = create_vessel_tube(x_min=0, x_max=VOLUME_WIDTH_IN_MM,
                                              z_min=0, z_max=VOLUME_HEIGHT_IN_MM,
                                              r_min=1, r_max=3)

    return tissue_dict

# Seed the numpy random configuration prior to creating the settings file in
# order to ensure that the same volume
# is generated with the same random seed every time.

np.random.seed(RANDOM_SEED)

settings = {
    # These parameters set the general properties of the simulated volume
    Tags.RANDOM_SEED: RANDOM_SEED,
    Tags.VOLUME_NAME: "MyVolumeName_" + str(RANDOM_SEED),
    Tags.SIMULATION_PATH: SAVE_PATH,
    Tags.SPACING_MM: SPACING,
    Tags.DIM_VOLUME_Z_MM: VOLUME_HEIGHT_IN_MM,
    Tags.DIM_VOLUME_X_MM: VOLUME_WIDTH_IN_MM,
    Tags.DIM_VOLUME_Y_MM: VOLUME_WIDTH_IN_MM,
    Tags.AIR_LAYER_HEIGHT_MM: 0,
    Tags.GELPAD_LAYER_HEIGHT_MM: 0,

    # The following parameters set the optical forward model
    Tags.RUN_OPTICAL_MODEL: True,
    Tags.WAVELENGTHS: np.arange(700, 951, 10),
    Tags.OPTICAL_MODEL_NUMBER_PHOTONS: 1e7,
    Tags.OPTICAL_MODEL_BINARY_PATH: MCX_BINARY_PATH,
    Tags.OPTICAL_MODEL: Tags.MODEL_MCX,
```

```

Tags.ILLUMINATION_TYPE: Tags.ILLUMINATION_TYPE_PENCIL,
Tags.LASER_PULSE_ENERGY_IN_MILLIJOULE: 50,

# The following parameters tell the script that we do not want any extra
# modelling steps
Tags.RUN_ACOUSTIC_MODEL: False,
Tags.APPLY_NOISE_MODEL: False,
Tags.PERFORM_IMAGE_RECONSTRUCTION: False,
Tags.SIMULATION_EXTRACT_FIELD_OF_VIEW: False,

# Add the structures to be simulated to the tissue
Tags.STRUCTURES: create_example_tissue()
}
print("Simulating ", RANDOM_SEED)
simulate(settings)
# TODO settings[Tags.IPPAI_OUTPUT_PATH]
print("Simulating ", RANDOM_SEED, "[Done]")

```

Reading the HDF5 simulation output

The file can be found in `samples/access_saved_PA1_data.py`:

```

from ippai.io_handling import load_hdf5, save_hdf5
import matplotlib.pyplot as plt
import numpy as np
from ippai.simulate import SaveFilePaths

PATH = "/home/janek/test/Vessels_10005/ippai_output.hdf5"
WAVELENGTH = 800 # currently only 800 and 900 are simulated as well

file = load_hdf5(PATH)

print(file['simulations'].keys())

fluence = (file['simulations']['original_data']['optical_forward_model_output']
           [str(WAVELENGTH)]['fluence'])
initial_pressure = (file['simulations']['original_data']
                    ['optical_forward_model_output']
                    [str(WAVELENGTH)]['initial_pressure'])
absorption = (file['simulations']['original_data']['simulation_properties']
              [str(WAVELENGTH)]['mua'])

shape = np.shape(fluence)

if len(shape) > 2:
    plt.figure()
    plt.subplot(231)
    plt.imshow(np.log10(fluence[int(shape[0]/2), :, :]))
    plt.subplot(232)
    plt.imshow(np.log10(absorption[int(shape[0]/2), :, :]))
    plt.subplot(233)
    plt.imshow(np.log10(initial_pressure))
    plt.subplot(234)
    plt.imshow(np.log10(fluence[:, int(shape[1]/2), :]))
    plt.subplot(235)
    plt.imshow(np.log10(absorption[:, int(shape[1]/2), :]))
    plt.subplot(236)
    plt.imshow(np.log10(initial_pressure))
    plt.show()
else:

```

```

plt.figure()
plt.subplot(131)
plt.imshow(np.log10(fluence[1:129, -65:-1]))
plt.subplot(132)
plt.imshow(np.log10(absorption[1:129, -65:-1]))
plt.subplot(133)
plt.imshow(np.log10(initial_pressure[1:129, -65:-1]))
plt.show()

save_hdf5(file, _PATH)

```

Defining custom tissue structures and properties

The file can be found in `samples/create_custom_tissues.py`:

```

from ippai.utils import TissueSettingsGenerator
from ippai.utils import CHROMOPHORE_LIBRARY
from ippai.utils import Chromophore
from ippai.utils import AbsorptionSpectrum
import numpy as np

def create_custom_absorber():
    wavelengths = np.linspace(200, 1500, 100)
    absorber = AbsorptionSpectrum(spectrum_name="random absorber",
                                wavelengths=wavelengths,
                                absorption_per_centimeter=np.random.random(
                                    np.shape(wavelengths)))

    return absorber

def create_custom_chromophore(volume_fraction: float = 1.0):
    chromophore = Chromophore(
        spectrum=create_custom_absorber(),
        volume_fraction=volume_fraction,
        musp500=40.0,
        b_mie=1.1,
        f_ray=0.9,
        anisotropy=0.9
    )
    return chromophore

def create_custom_tissue_type():
    # First create an instance of a TissueSettingsGenerator
    tissue_settings_generator = TissueSettingsGenerator()

    water_volume_fraction = 0.4
    bvf = 0.5
    oxy = 0.4

    # Then append chromophores that you want
    tissue_settings_generator.append(key="oxyhemoglobin", value=
                                    CHROMOPHORE_LIBRARY.oxyhemoglobin(oxy*bvf))
    tissue_settings_generator.append(key="deoxyhemoglobin", value=
                                    CHROMOPHORE_LIBRARY.deoxyhemoglobin(oxy * bvf))
    tissue_settings_generator.append(key="water", value=
                                    CHROMOPHORE_LIBRARY.water(water_volume_fraction))
    tissue_settings_generator.append(key="custom", value=

```



```

        create_custom_chromophore(0.1))

return tissue_settings_generator.get_settings()

```

Class references

Module: utils

class `ippai.utils.libraries.literature_values.MorphologicalTissueProperties`

This class contains a listing of morphological tissue parameters as reported in literature. The listing is not the result of a meta analysis, but rather uses the best fitting paper at the time of implementation. Each of the fields is annotated with a literature reference or a descriptions of how the particular values were derived for tissue modelling.

class `ippai.utils.libraries.literature_values.OpticalTissueProperties`

This class contains a listing of optical tissue parameters as reported in literature. The listing is not the result of a meta analysis, but rather uses the best fitting paper at the time of implementation. Each of the fields is annotated with a literature reference or a descriptions of how the particular values were derived for tissue modelling.

class `ippai.utils.libraries.literature_values.StandardProperties`

This class contains a listing of default parameters that can be used. These values are sensible default values but are generally not backed up by proper scientific references, or are rather specific for internal use cases.

class `ippai.utils.libraries.spectra_library.AbsorptionSpectrum` (`spectrum_name: str`,
`wavelengths: numpy.ndarray`, `absorption_per_centimeter: numpy.ndarray`)

An instance of this class represents the absorption spectrum over wavelength for a particular

`get_absorption_for_wavelength(wavelength: float) → float`

Parameters: `wavelength` – the wavelength to retrieve a optical absorption value for [cm^{-1}].

Returns: the best matching linearly interpolated absorption value for the given wavelength.

`get_absorption_over_wavelength()`

Returns: numpy array with the available wavelengths and the corresponding absorption properties

`ippai.utils.libraries.spectra_library.view_absorption_spectra` (`save_path=None`)

Opens a matplotlib plot and visualizes the available absorption spectra.

Parameters: `save_path` – If not None, then the figure will be saved as a png file to the destination.

class `ippai.utils.libraries.tissue_library.TissueLibrary`

TODO

`blood_arterial()`

Returns: a settings dictionary containing all min and max parameters fitting for full blood.

`blood_generic()`

Returns: a settings dictionary containing all min and max parameters fitting for full blood.

`blood_venous()`

Returns: a settings dictionary containing all min and max parameters fitting for full blood.

`bone()`

Returns: a settings dictionary containing all min and max parameters fitting for full blood.

```
constant (mua, mus, g)
    TODO
```

```
dermis (background_oxy=0.5)
```

Returns: a settings dictionary containing all min and max parameters fitting for dermis tissue.

```
epidermis ()
```

Returns: a settings dictionary containing all min and max parameters fitting for epidermis tissue.

```
get_blood_volume_fractions (total_blood_volume_fraction, oxygenation)
    TODO
```

```
muscle (background_oxy=0.5)
```

Returns: a settings dictionary containing all min and max parameters fitting for generic background tissue.

```
subcutaneous_fat (background_oxy=0.5)
```

Returns: a settings dictionary containing all min and max parameters fitting for subcutaneous fat tissue.

```
class ippai.utils.libraries.tissue_library.TissueSettingsGenerator
    TODO
```

Module: io_handling

```
ippai.io_handling.io_hdf5.load_hdf5 (file_path, file_dictionary_path='/')
```

Loads a dictionary from an hdf5 file. :param file_path: Path of the file to load the dictionary from. :param file_dictionary_path: Path in dictionary structure of hdf5 file to load the dictionary in. :return: Dictionary

```
ippai.io_handling.io_hdf5.save_hdf5 (dictionary, file_path, file_dictionary_path='/')
```

Saves a dictionary with arbitrary content to an hdf5-file with given filepath. :param dictionary: Dictionary to save. :param file_path: Path of the file to save the dictionary in. :param file_dictionary_path: Path in dictionary structure of existing hdf5 file to store the dictionary in. return

Module: simulate

```
ippai.simulate.simulation.simulate (settings)
```

Parameters: settings –

Returns:

```
class ippai.simulate.tissue_properties.TissueProperties (settings: dict)
```

```
ensure_valid_settings (settings: dict)
```

Parameters: settings – a dictionary containing at least one key - value pair of a string and a corresponding Chromophore instance.

```
get (wavelength)
    TODO
```

```
ippai.simulate.volume_creator.create_acoustic_properties (volumes, settings)
```

Creates maps of density, speed of sound and acoustic attenuation based on the segmented mask in volumes.

Parameters:

- **volumes** – The volumes to append the acoustic parameters to
- **settings** – The settings to extract if the medium is homogeneous or if heterogeneous parameters should be set.

Returns: The volumes with appended maps of acoustic parameters of size `volumes[0].size`

```
ippai.simulate.volume_creator.create_gruneisen_map(volumes, settings)
```

Creates a map the gruneisenparameter based on the temperature given in `Tags.MEDIUM_TEMPERATURE_CELCIUS`. If no medium temperature is specified, then a standard body temperature of 36°C is assumed.

Parameters:

- **volumes** – The volumes to append the gruneisen parameter to
- **settings** – The settings to extract the temperature from

Returns: the volumes with an appended map of gruneisen parameters of size `volumes[0].size`

```
ippai.simulate.volume_creator.create_simulation_volume(settings)
```

This method creates a in silico representation of a tissue as described in the settings file that is given. :param settings: a dictionary containing all relevant Tags for the simulation to be able to instantiate a tissue. :return: a path to a npz file containing characteristics of the simulated volume:

absorption, scattering, anisotropy, oxygenation, and a segmentation mask. All of these are given as 3d numpy arrays.

```
ippai.simulate.volume_creator.fnc_straight_ellipse(x, y, z, r_x, r_z, X1, X2)
```

cartesian representation of a straight tube that goes from position X1 to position X2 with radius r. :param x: :param y: :param z: :param r: :param X1: :param X2: :return:

```
ippai.simulate.volume_creator.fnc_straight_tube(x, y, z, r, X1, X2)
```

cartesian representation of a straight tube that goes from position X1 to position X2 with radius r. :param x: :param y: :param z: :param r: :param X1: :param X2: :return:

```
ippai.simulate.volume_creator.merge_voxel(volumes, x_idx, y_idx, z_idx, mua, mus, g, oxy, seg, fraction)
```

Updates a voxel position in the volumes by merging the given physical properties with the properties already stored in the volumes. The merging is done in a relative manner using the given fraction.

Parameters:

- **volumes** – list of numpy arrays with `len(volumes) >= 3`
- **x_idx** – integer
- **y_idx** – integer
- **z_idx** – integer
- **mua** – scalar, the optical absorption coefficient in 1/cm
- **mus** – scalar, the optical scattering coefficient in 1/cm
- **g** – scalar, the anisotropy
- **oxy** – scalar, the blood oxygenation in [0, 1]
- **seg** – integer, the tissue segmentation type from `SegmentationClasses`
- **fraction** – scalar in [0, 1]

Returns: the volumes with the changed properties

```
ippai.simulate.volume_creator.set_voxel(volumes, x_idx, y_idx, z_idx, mua, mus, g, oxy, seg)
```

Sets a voxel position to a specific value in the volume

Parameters:

- **volumes** – list of numpy arrays with `len(volumes) >= 3`
- **y_idx** – integer
- **x_idx** – integer
- **z_idx** – integer
- **mua** – scalar
- **mus** – scalar
- **g** – scalar

Returns: the volumes with the changed properties

Conventions

`class ippai.simulate.constants.SaveFilePaths`

The save file paths specify the path of a specific data structure in the dictionary of the ippai output hdf5. All of these paths have to be used like: `SaveFilePaths.PATH.format(Tags.UPSAMPLED_DATA or Tags.ORIGINAL_DATA, wavelength)`

`class ippai.simulate.constants.SegmentationClasses`

The segmentation classes define which “tissue types” are modelled in the simulation volumes.

Structures

`ippai.simulate.structures.create_subcutaneous_fat_layer(background_oxy=0.0)`

FIXME: DEPRECATED FIXME: first research into fat tissue properties before adding this to simulation (!) :param background_oxy: :return:

Forward models

`class ippai.simulate.models.optical_models.OpticalForwardAdapterBase`

Use this class as a base for implementations of optical forward models.

abstract forward_model (absorption_cm, scattering_cm, anisotropy, settings)

A deriving class needs to implement this method according to its model.

Parameters:

- **absorption_cm** – Absorption in units of per centimeter
- **scattering_cm** – Scattering in units of per centimeter
- **anisotropy** – Dimensionless scattering anisotropy
- **settings** – Setting dictionary

Returns: Fluence in units of J/cm²

simulate (optical_properties_path, settings)

Call this method to invoke the simulation process.

A adapter that implements the forward_model method, will take optical properties of absorption, scattering, and scattering anisotropy as input and return the light fluence as output.

Parameters:

- **optical_properties_path** – path to a *.npz file that contains the following tags:
Tags.PROPERTY_ABSORPTION_PER_CM -> contains the optical absorptions in units of one per centimeter
Tags.PROPERTY_SCATTERING_PER_CM -> contains the optical scattering in units of one per centimeter
Tags.PROPERTY_ANISOTROPY -> contains the dimensionless optical scattering anisotropy
- **settings** –

Returns:

`class ippai.simulate.models.acoustic_models.AcousticForwardAdapterBase`
Use this class as a base for implementations of optical forward models.

`abstract forward_model(settings)`

A deriving class needs to implement this method according to its model.

Parameters: **settings** – Setting dictionary

Returns: Fluence in units of J/cm²

`simulate(optical_properties_path, settings)`

Call this method to invoke the simulation process. TODO

A adapter that implements the forward_model method, will take acoustic properties as input and return the time series pressure data as output.

Parameters:

- **optical_properties_path** – path to a *.npz file that contains the following tags:
Tags.PROPERTY_ABSORPTION_PER_CM -> contains the optical absorptions in units of one per centimeter
Tags.PROPERTY_SCATTERING_PER_CM -> contains the optical scattering in units of one per centimeter
Tags.PROPERTY_ANISOTROPY -> contains the dimensionless optical scattering anisotropy

- **settings** –

Returns:

`class ippai.simulate.models.noise_models.GaussianNoiseModel`

This class is responsible to apply an additive gaussian noise to the input data.

`apply_noise_model(data, settings)`

Parameters:

- **data** –
- **settings** –

Returns:

Index

A

[AbsorptionSpectrum](#) (class in [ippai.utils.libraries.spectra_library](#))

[AcousticForwardAdapterBase](#) (class in [ippai.simulate.models.acoustic_models](#))

[apply_noise_model\(\)](#) ([ippai.simulate.models.noise_models.GaussianNoiseModel](#) method)

B

[blood_arterial\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

[blood_generic\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

[blood_venous\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

[bone\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

C

[constant\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

[create_acoustic_properties\(\)](#) (in module [ippai.simulate.volume_creator](#))

[create_gruneisen_map\(\)](#) (in module [ippai.simulate.volume_creator](#))

[create_simulation_volume\(\)](#) (in module [ippai.simulate.volume_creator](#))

[create_subcutaneous_fat_layer\(\)](#) (in module [ippai.simulate.structures](#))

D

[dermis\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

E

[ensure_valid_settings\(\)](#) ([ippai.simulate.tissue_properties.TissueProperties](#) method)

[epidermis\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

F

[fnc_straight_ellipse\(\)](#) (in module [ippai.simulate.volume_creator](#))

[fnc_straight_tube\(\)](#) (in module [ippai.simulate.volume_creator](#))

[forward_model\(\)](#) ([ippai.simulate.models.acoustic_models.AcousticForwardAdapterBase](#) method)

[forward_model\(\)](#) ([ippai.simulate.models.optical_models.OpticalForwardAdapterBase](#) method)

G

[GaussianNoiseModel](#) (class in [ippai.simulate.models.noise_models](#))

[get\(\)](#) ([ippai.simulate.tissue_properties.TissueProperties](#) method)

[get_absorption_for_wavelength\(\)](#) ([ippai.utils.libraries.spectra_library.AbsorptionSpectrum](#) method)

[get_absorption_over_wavelength\(\)](#) ([ippai.utils.libraries.spectra_library.AbsorptionSpectrum](#) method)

[get_blood_volume_fractions\(\)](#) ([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

I

[ippai.io_handling](#) (module)

[ippai.io_handling.io_hdf5](#) (module)

[ippai.simulate](#) (module)

[ippai.simulate.constants](#) (module)

[ippai.simulate.models](#) (module)

[ippai.simulate.models.acoustic_models](#) (module)

[ippai.simulate.models.noise_models](#) (module)

[ippai.simulate.models.optical_models](#) (module)

[ippai.simulate.models.reconstruction_models](#) (module)

[ippai.simulate.simulation](#) (module)

[ippai.simulate.structures](#) (module)

[ippai.simulate.tissue_properties](#) (module)

[ippai.simulate.volume_creator](#) (module)

[ippai.utils](#) (module)

[ippai.utils.libraries](#) (module)

[ippai.utils.libraries.chromophore_library](#) (module)

[ippai.utils.libraries.literature_values](#) (module)

[ippai.utils.libraries.spectra_library](#) (module)

[ippai.utils.libraries.tissue_library](#) (module)

L

[load_hdf5\(\)](#) (in module [ippai.io_handling.io_hdf5](#))

M

[merge_voxel\(\)](#) (in module [ippai.simulate.volume_creator](#))

[MorphologicalTissueProperties](#) (class in [ippai.utils.libraries.literature_values](#))

[muscle\(\)](#)
([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

O

[OpticalForwardAdapterBase](#) (class in [ippai.simulate.models.optical_models](#))

[OpticalTissueProperties](#) (class in [ippai.utils.libraries.literature_values](#))

S

[save_hdf5\(\)](#) (in module [ippai.io_handling.io_hdf5](#))

[SaveFilePaths](#) (class in [ippai.simulate.constants](#))

[SegmentationClasses](#) (class in [ippai.simulate.constants](#))

[set_voxel\(\)](#) (in module [ippai.simulate.volume_creator](#))

[simulate\(\)](#) (in module [ippai.simulate.simulation](#))

([ippai.simulate.models.acoustic_models.AcousticForwardAdapterBase](#) method)

([ippai.simulate.models.optical_models.OpticalForwardAdapterBase](#) method)

[StandardProperties](#) (class in [ippai.utils.libraries.literature_values](#))

[subcutaneous_fat\(\)](#)
([ippai.utils.libraries.tissue_library.TissueLibrary](#) method)

T

[TissueLibrary](#) (class in [ippai.utils.libraries.tissue_library](#))

[TissueProperties](#) (class in [ippai.simulate.tissue_properties](#))

[TissueSettingsGenerator](#) (class in [ippai.utils.libraries.tissue_library](#))

V

[view_absorption_spectra\(\)](#) (in module [ippai.utils.libraries.spectra_library](#))

Python Module Index

i

[ippai](#)

[ippai.io_handling](#)

[ippai.io_handling.io_hdf5](#)

[ippai.simulate](#)

[ippai.simulate.constants](#)

[ippai.simulate.models](#)

[ippai.simulate.models.acoustic_models](#)

[ippai.simulate.models.noise_models](#)

[ippai.simulate.models.optical_models](#)

[ippai.simulate.models.reconstruction_models](#)

[ippai.simulate.simulation](#)

[ippai.simulate.structures](#)

[ippai.simulate.tissue_properties](#)

[ippai.simulate.volume_creator](#)

[ippai.utils](#)

[ippai.utils.libraries](#)

[ippai.utils.libraries.chromophore_library](#)

[ippai.utils.libraries.literature_values](#)

[ippai.utils.libraries.spectra_library](#)

[ippai.utils.libraries.tissue_library](#)